6 X 0,7 ML PREFILLED SYRINGES

# view ILM

0,03% BLULIFE

view ILM is 0.03% blulife

## **TISSUE STAINING**



	view ILM
Ingredients	<ul> <li>blulife (C<sub>48</sub>H<sub>50</sub>N<sub>3</sub>NaO<sub>7</sub>S)</li> <li>thickening agent</li> <li>sodium chloride</li> <li>sodium phosphate monobasic</li> <li>sodium phosphate bibasic</li> <li>water for injectables</li> </ul>
Dye concentration	0,027 -0,033% w/v
Hydrogenionic activity (pH)	7,00 – 7,40 pH units
Osmolality	270 – 330 mOsm/kg
Bacterial endotoxins	< 0,5 EU/ml
Form	liquid
Volume	≥ 0,7 ml
Packaging	syringe
Sterilization	moist heat according to EN 556-1 and EN ISO 17665-1
Storage temperature	+15°C / +32°C
Expiry date	36 months
Intended use	reversible staining of the inner limiting membrane (ILM) during ophthalmic surgery
Class	medical device
GMDN code	P 45180
CND code	Q 02030201
MD registration number	1331840/R
Medical device class	lla
CE mark	CE 0546

 $(C_{48}H_{50}N_3NaO_7S)$ physiological solution, prefilled packaged in syringes of 0,7 ml each, for reversible staining of the inner limiting membrane (ILM). The product warrants a prompt and effective tissue staining, without any risks for the patient. It contains a thickening agent that improves the contact between the dye and the tissue, thus allowing a clear visualization of the inner limiting membrane. In addition, it simplifies the surgery, making superfluous the fluid/ air exchange performed before the dye injection, in order not to dilute it.

#### Ref.: RMB 003-00

#### **Contents:**

- ▶ 6 view ILM prefilled syringes, 0,7 ml
- ▶ 6 pistons
- ▶ 1 insert
- ▶ 12 labels for product traceability

GS1 – 128 (EAN 128) barcode including EAN13 barcode:

8021118001132

### **ALCHIMIA**

Viale Austria, 14 35020 Ponte S.Nicolò (PD) Italy

Tel. +39 049 8962074 / 64 Fax +39 049 8962071

www.alchimiasrl.com info@alchimiasrl.com