

MEDIA FOR  
DONOR GLOBE  
AND CORNEAL  
RINSING

# PSS-L

PHOSPHATE BUFFERED SALINE SOLUTION

## GLOBE AND CORNEAL RINSING



PSS-L is a phosphate buffered saline solution for donor globe rinsing, before corneal excision, and for corneal rinsing during trypan blue evaluation.

PSS-L	
<b>Ingredients</b>	<ul style="list-style-type: none"> <li>▶ sodium chloride</li> <li>▶ sodium phosphate</li> <li>▶ potassium chloride</li> <li>▶ potassium phosphate</li> <li>▶ sodium hydroxide</li> </ul>
<b>Form</b>	sterile liquid
<b>Volume</b>	30 ml
<b>Physico-chemical parameters</b>	<ul style="list-style-type: none"> <li>▶ hydrogenionic activity (pH) : 7.40 (7.20 – 7.60)</li> <li>▶ chloride content (% p/p) : 0.80 – 0.88</li> <li>▶ phosphate identification : positive</li> <li>▶ osmolality (MOsm/Kg) : 300 (255 - 345)</li> </ul>
<b>Sterilization</b>	sterile (filtration and aseptic distribution according to UNI EN 556-2 and UNI EN ISO 13408-1)
<b>Storage temperature</b>	+2°C / +25°C
<b>Expiry date</b>	48 months
<b>Intended use</b>	donor globe and corneal rinsing
<b>Temperature of actual use</b>	room temperature
<b>Class</b>	medical device
<b>GMDN code</b>	P36094
<b>MD registration number</b>	3736/R
<b>Medical device class</b>	Ila
<b>CE mark</b>	CE 0546

**Ref.:** GRS 003-00

**Contents:**

- ▶ 12 30ml-vial pack

**EAN 13 barcode:**

8021118000463

## ALCHIMIA

Viale Austria, 14  
35020 Ponte S.Nicolò (PD)  
Italy

P.IVA IT 00063370282

Phone +39 049 8962074 / 64  
Fax +39 049 8962071

www.alchimiasrl.com  
info@alchimiasrl.com