

**6 X 0.7 ML PREFILLED SYRINGES** 

# **TISSUE STAINING**

	viewILM
Ingredients	<ul> <li>blulife (C<sub>48</sub>H<sub>50</sub>N<sub>3</sub>NaO<sub>7</sub>S)</li> <li>thickening agent</li> <li>sodium chloride</li> <li>sodium phosphate monobasic</li> <li>sodium phosphate bibasic</li> <li>water for injectables</li> </ul>
Dye concentration	0.027 - 0.033 %w/V
Hydrogenionic activity (pH)	7.00 - 7.40 pH units
Osmolality	245 - 360 mOsm/kg
Density	1.015 - 1.018 g/ml
Bacterial Endotoxins	< 0.5 EU/ml
Cytotoxicity Test	Non-cytotoxic according to ISO 10993-5, direct contact test performed on two cell lines: - retinal pigment epithelial cells (ARPE-19) - mouse fibroblasts cells (BALB 3T3 clone A31)
Form	sterile liquid
Volume	≥ 0.7 ml
Packaging	syringe
Sterilization	moist heat according to EN 556-1 and EN ISO 17665-1
Storage temperature	+15°C / +32°C
Expiry date	36 months
Intended use	reversible staining of the inner limiting membrane (ILM) during ophthalmic surgery
Class	medical device
GMDN code	P 45180
CND code	Q 02030201
MD registration number	1331840/R
Medical device class	lla
CE mark	CE 0477

viewILM is 0.03% blulife (C48H50N3NaOzS) in physiological solution, packaged in prefilled syringes of 0.7 ml each, for reversible staining of the inner limiting membrane (ILM). The product warrants a prompt and effective tissue staining, without any risks for the patient. It contains a thickening agent that improves the contact between the dye and the tissue, thus allowing a clear visualization of the inner limiting membrane. In addition, it simplifies the surgery, making superfluous the fluid/air exchange performed before the dye injection, in order not to dilute it.

### Ref.: RMB 003-00

#### **Contents:**

- 6 viewILM prefilled syringes, 0.7 ml
- ▶ 6 pistons
- I package insert
- 12 labels for product traceability

#### GS1-128 barcode

# ALCHIMIA

Viale Austria, 14 35020 Ponte S.Nicolò (PD) Italy

VAT IT 00063370282

Phone +39 049 8962074 / 64 Fax +39 049 8962071

www.alchimiasrl.com info@alchimiasrl.com