

6 X 0,7 ML
PREFILLED SYRINGES

TWIN

0,18% TRYPAN BLUE + 0,03% BLULIFE

TISSUE STAINING



TWIN	
Ingredienti	<ul style="list-style-type: none"> ▶ trypan blue ▶ blulife (C₄₈H₅₀N₃NaO₇S) ▶ thickening agent ▶ sodium chloride ▶ sodium phosphate monobasic ▶ sodium phosphate bibasic ▶ water for injectables
Dye concentration	0,162 – 0,198% w/v trypan blue 0,027 – 0,033% w/v C ₄₈ H ₅₀ N ₃ NaO ₇ S
Hydrogenionic activity (pH)	7,00 – 7,40 pH units
Osmolality	270 – 330 mOsm/kg
Bacterial endotoxins	< 0,5 EU/ml
Form	liquid
Volume	≥ 0,7 ml
Packaging	syringe
Sterilization	moist heat according to EN 556-1 and EN ISO 17665-1
Storage temperature	+15°C / +32°C
Expiry date	36 months
Intended use	reversible staining of the retinal membranes and the inner limiting membrane (ILM) during ophthalmic surgery
Class	medical device
GMDN code	P 45180
CND code	Q 02030201
MD registration number	1331847/R
Medical device class	Ila
CE mark	CE 0546

TWIN is 0,18% trypan blue + 0,03% blulife (C₄₈H₅₀N₃NaO₇S) in physiological solution, packaged in prefilled syringes of 0,7 ml each, for reversible staining of the retinal membranes and the inner limiting membrane (ILM). The product warrants a prompt and effective tissue staining, without any risks for the patient. It contains a thickening agent that improves the contact between the dyes and the tissues, thus allowing a clear visualization of the retinal membranes and the inner limiting membrane. In addition, it simplifies the surgery, making superfluous the fluid/air exchange performed before the dye injection, in order not to dilute it.

Ref.: RMB 004-00

Contents:

- ▶ 6 TWIN prefilled syringes, 0,7 ml
- ▶ 6 pistons
- ▶ 1 insert
- ▶ 12 labels for product traceability

**GS1 – 128 (EAN 128)
barcode including
EAN13 barcode:**

8021118001149

ALCHIMIA

Viale Austria, 14
35020 Ponte S.Nicolò (PD)
Italy

Tel. +39 049 8962074 / 64
Fax +39 049 8962071

www.alchimiasrl.com
info@alchimiasrl.com