

# OCIGEL

10 SINGLE-DOSE 5 UNITS STRIPS OF 3 ML

## COUPLING GEL FOR OCULAR DIAGNOSTICS BY CONTACT

	OCIGEL
<b>Ingredients</b>	<ul style="list-style-type: none"> <li>▶ carbopol® 980</li> <li>▶ propylene glycol</li> <li>▶ purified water</li> </ul>
<b>Form</b>	sterile gel
<b>Volume</b>	3 ml
<b>Chemical-Physical Parameters</b>	
<b>Hydrogenionic activity</b>	7.00 - 7.50 pH units
<b>Osmolality</b>	265 - 320 mOsm/kg
<b>Dynamic Viscosity</b>	66 - 100 Pa*s
<b>Transmittance</b>	≥ 90%
<b>Permanence on lens</b>	> 30 min
<b>Cytotoxicity test</b>	non-cytotoxic in accordance with EN ISO 10993-5 (direct contact test, performed on BALB 3T3 cell line)
<b>Packaging</b>	single-dose unit
<b>Sterilization</b>	steam sterilization according to UNI EN 556-1 and EN ISO 17665-1
<b>Storage temperature</b>	room temperature
<b>Expiry date</b>	30 months
<b>Intended use</b>	coupling gel for ocular diagnostics by contact
<b>Category</b>	medical device
<b>CND code</b>	Q02030399
<b>Device class</b>	I sterile
<b>CE mark</b>	CE 0546

Viscous and transparent gel, non-cytotoxic, preservative-free, with high ocular tolerability, intended to be used as an interface (coupling gel) with ocular diagnostics instruments by contact; for use in hospitals and medical offices.

Simple to apply with no air bubbles and easy to remove after use.

Offers high stability on contact lenses for diagnostics and laser therapy.

Available in packs of 50 sterile single-dose units of 3 ml each.

**Ref.: OCI 002**

**Contents:**

- ▶ 10 strips of 5 sterile single-dose units of 3 ml
- ▶ 1 package insert

**GS1-128 barcode**

**ALCHIMIA**

Viale Austria, 14  
35020 Ponte S.Nicolò (PD)  
Italy

VAT IT 00063370282

Phone +39 049 8962074 / 64  
Fax +39 049 8962071

www.alchimiasrl.com  
info@alchimiasrl.com